

Pyrophosphatase, Inorganic, GMP-Grade

Catalog #GMP-PYR-YE101

Storage Condition -20°C ± 5°C for 24 months. Avoid repeated freeze/thaw cycles.

Form Liquid

Source An *E. coli* strain that carries the gene for pyrophosphatase, Inorganic from *Saccharomyces cerevisiae*

Storage Buffer 20 mM Tris-HCl, 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, and 50% Glycerol, pH 8.0

Concentration 0.1U/μL

Unit Definition One unit is defined as the amount of enzyme that will generate 1μmol of orthophosphate per minute from inorganic pyrophosphate (PPi) in a 10 min standard reaction at 25°C containing 20 mM Tris-HCl, pH 8.0, 2 mM MgCl₂, and 2 mM PPi.

Product Contents

- Pyrophosphatase, Inorganic (0.1U/μL)

Product Description

Pyrophosphatase, Inorganic catalyzes the hydrolysis of inorganic pyrophosphate to form two orthophosphates. Pyrophosphatase, Inorganic enhances RNA or DNA synthesis by preventing the accumulation of byproduct PPi generated in the reaction.

Applications

- Improving RNA yield in *In vitro* transcription (IVT)
- Enhancing DNA amplification reaction
-

Quality Control Statement

KACTUS manufactures this product according to GMP guidelines and performs stringent quality control testing before release. It is suitable for mRNA therapeutics and manufacturing of mRNA vaccines. Regulatory support documents are available. Please contact support@kactusbio.us for more information.

Quality Control Release Criteria

Assay	Criteria
Activity (Enzyme Catalytic Efficiency)	≥ 126U/mL
Purity (SEC-HPLC)	≥ 95%
Residual Endonuclease	Negative
Residual Exonuclease	Negative
Endotoxin	≤ 10EU/mL
Residual RNase	Negative
Residual DNase	Negative
Residual Protease	Negative

Residual Host Cell Protein	≤ 20ng/mg
Residual Host Cell DNA	≤ 100pg/mL
Bioburden	≤ 1CFU/10mL

Protocol for *In Vitro* Transcription

1. Prepare the following reaction mixture:


Reagent	Quantity
5X Transcription Buffer-1 (included with T7 RNA Polymerase)	4μL
CTP/GTP/ATP/UTP (100mM each)	2μL each
Murine RNase Inhibitor	1μL
Pyrophosphatase, Inorganic	1μL
T7 RNA Polymerase	2μL
Template DNA	1μg
RNase-free Water	Up to 20μL

2. Incubate at 37°C for 1-2 hours.
3. After transcription, add 2U [DNase I](#) to digest DNA template for 15 minutes at 37°C.
4. Inactivate DNase I by phenol/chloroform extraction.

Notes

- Pyrophosphatase, Inorganic works well in any buffer containing Mg²⁺ (1-10 mM)
- The hydrolysis of PPi is a proportional to the concentration of enzyme.

Performance Validation



The addition of Pyrophosphatase, Inorganic to IVT reactions improves mRNA yield. Lane 1: Without Pyrophosphatase, Inorganic; Lane 2: With Pyrophosphatase, Inorganic.