

Mycoplasma DNA (qPCR) Detection Kit

CATALOG #DK-M002, 50 Tests

PRODUCT INTRODUCTION

The Mycoplasma DNA (qPCR) Detection Kit is designed for qualitative detection of Mycoplasma contamination in critical biological materials, including master cell banks, working cell banks, viral seed lots, control cells, and cells for clinical treatment. Utilizing TaqMan® probe-based qPCR technology, this kit enables quantitative detection of Mycoplasma DNA in samples, delivering rapid, specific, and reliable results while covering over 150 Mycoplasma species.

It is recommended to use KACTUS Host Cell Residual DNA (Magnetic Beads) Sample Preparation Kit (Catalog No. DK-SP001) for integrated workflow solutions.

KIT COMPONENTS

Components	Size	Storage Conditions
Mycoplasma Primer & Probe MIX	120 µL × 1 vial	-20°C, protected from light
2X qPCR Reaction MIX	700 µL × 1 vial	-20°C
Internal Control (IC)	220 µL × 1 vial	
Positive Control (PC)	500 µL × 1 vial	
DNA Diluent	1.5 mL × 3 vials	

STORAGE AND SHELF LIFE

Unopened reagents maintain stability for 12 months from the date of manufacture when stored as recommended. After opening, please store the remaining reagents at -20°C.

REQUIRED MATERIALS NOT SUPPLIED

- Real-time PCR instrument (ABI QuantStudio 3, ABI 7500, etc.), Centrifuge, Vortex.
- Optical 8-tube strip or 96-well qPCR plate (sterile, nuclease-free, and low retention), 1.5mL centrifuge tube (sterile and low retention)

WORKFLOW

Prepare the reagents and samples

1. **Prepare the kit.** Place all components on ice to ensure complete dissolution before use.
2. **Extract DNA from test samples.** It is recommended to extract sample DNA using KACTUS Host Cell Residual DNA (Magnetic Beads) Sample Pretreatment Kit (Catalog No. DK-SP001).
3. **Prepare the qPCR Working Mix.**
 - a) Calculate the required number of qPCR reaction wells using the following formula:

- i. Total wells N = (1 Non-Template Control + 1 Positive Control + Test Samples) × 2 or 3 Number of Replicates
- b) Prepare the Working Mix per table below:
 - i. Total volume of Working Mix = 15 µL × (Total Wells N + 2)

Components	Working Mix for single reaction
2×qPCR Reaction MIX	12.5 µL
Mycoplasma Primer & Probe MIX	2.1 µL
Internal Control (IC)	0.4 µL
Total volume	15 µL

- c) Prepare the qPCR reaction solution. Place all the reagents on ice, prepare the qPCR reaction solution per the table below:


	Sample	Positive Control (PC)	Non-Template Control (NTC)
Template	10 µL Sample	10 µL Positive Control (PC)	10 µL DNA Dilution
Working Mix per well	15 µL	15 µL	15 µL
Total Volume per well	25 µL	25 µL	25 µL

qPCR reaction parameters

Using American ABI 7500 Real-Time PCR System (version 2.06) as an example:

1. Create a blank new program. Experiment Type: Quantitation-Standard Curve. Reagents: TaqMan® Reagents.
2. Use “Add New Targets” to create new probes and edit the probe names in “Target Name”. Select the desired fluorescent and quenching groups in “Reporter” and “Quencher”. “Reporter” is selected as FAM, the “Quencher” is selected as NONE, the fluorescent group of internal control (IC) is selected as VIC, and the reference fluorescence is ROX.
3. Use “Add New Samples” to add new samples and name them in the corresponding Sample name column.
4. Set the reaction parameters. 58°C for 34 s in Stage 3 is set as fluorescence collection.

Stage 1	Digest	Reps: 1	37°C	2 min
Stage 2	Pre-denaturation	Reps: 1	95°C	30 s
Stage 3	Cycle	Reps: 45	95°C	10 s
			58°C	34 s



5. Click “Start Run” to perform qPCR.

ANALYSIS

Using American ABI 7500 Real-Time PCR System (version 2.06) as an example:

1. Click “Analyze”. Check whether the shape of the amplification curve is normal.
2. Read the Ct value in the “View Well Table”.
3. NTC, PC results must meet:

Group	FAM channel	VIC channel
NTC	Ct ≥ 40 or no significant amplification curve observed.	Ct < 40 with a characteristic sigmoidal amplification curve.
PC	Ct < 40 with a characteristic sigmoidal amplification curve.	Ct < 40 with a characteristic sigmoidal amplification curve.

4. Test sample results:

Ct < 40 with a characteristic sigmoidal amplification curve.	Ct < 40 with a characteristic sigmoidal amplification curve.	Positive
	Ct ≥ 40 or no significant amplification curve observed.	PCR inhibition, retest.
Ct ≥ 40 or no significant amplification curve	Ct < 40 with a characteristic sigmoidal amplification curve.	Negative
	Ct ≥ 40 or no significant amplification curve	PCR inhibition, retest.

NOTES

1. The components in the kit are stable and do not require aliquoting.
2. Prepare negative and positive samples in separate areas avoiding cross-contamination.
3. Always use clean pipette tips and avoid contamination.
4. Use the kit before expiration and do not mix different lots.
5. Thaw all components completely before use. Mix well and centrifuge briefly.
6. Follow instructions carefully for the best results.
7. Results depend on reagent quality, technique, and lab conditions.
8. Prepare enough samples for your needs.
9. For research use only. Not for medical diagnosis.

FAQs

Common Issues	Possible Reasons	Actions to Take
No Ct value	Incorrect fluorescence detection settings in the PCR program.	Ensure 58°C for 34 s in Stage 3 is set as fluorescence collection.
	Primer or probe degradation	Perform PAGE electrophoresis to determine if degradation had occurred.
	Template insufficient.	Template degradation may be caused by impurities or damaged by repeated freeze-thaw.
Delayed Ct	qPCR reaction components degraded or are insufficient	Verify the ROX reference signal.
	Inhibitors in template	
Negative control shows a signal	Contamination	<ol style="list-style-type: none"> 1. Repeat using newly prepared Working Mix. 2. Prepare the reaction solution in the laminar flow hood. 3. Separate the different experiment areas. Change to low-retention, nuclease-free tubes and pipette tips.
Abnormal amplification curve	<ol style="list-style-type: none"> 1. High concentration of template; template degraded 2. Evaporation or sealing issue 3. Bubbles in tube 	<ol style="list-style-type: none"> 1. Dilute template. Retest with fresh reagents 2. Carefully transfer the reagents to avoid bubbles.

RELATED PRODUCTS

Product Name	Catalog No.
Host Cell Residual DNA (magnetic beads) Sample Preparation Kit	DK-SP001
Mycoplasma DNA (qPCR) Detection Kit	DK-M002