Human Activin A Protein

Cat. No. ACV-HM001




Cat. No. ACV-III	11001
Description	
Source	Recombinant Human Activin A Protein is expressed from HEK293 without tag.
	It contains Gly311-Ser426.
Accession	P08476
Molecular Weight	The protein has a predicted MW of 12.97 kDa. Due to glycosylation, the protein migrates to 14-15 kDa under reduced (R) condition, and 23-25 kDa under Non reducing (N) condition based on Tris-Bis PAGE result.
Endotoxin	Less than 0.1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in 4mM HCI. Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 μ g/ml is recommended. Dissolve the lyophilized protein in 4mM HCI.
Storage	-20 to -80°C for 12 months as supplied from date of receipt80°C for 3-6 months after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Activin-A, a member of the TGF-β superfamily, Extensive research over the past decades illuminated fundamental roles for activin-A in essential biologic processes, including embryonic development, stem cell maintenance and differentiation, haematopoiesis, cell proliferation and tissue fibrosis. Activin-A signals through two type I and two type II receptors which, upon ligand binding, activate their kinase activity, phosphorylate the SMAD2 and 3 intracellular signaling mediators that form a complex with SMAD4, translocate to the nucleus and activate or silence gene expression.

Assay Data

Tris-Bis PAGE



ELISA Data

Human Activin A on Tris-Bis PAGE under reduced (R) condition and Non reducing (N) condition. The purity is greater than 95%.

KAGTUS

Assay Data

Human Activin A, No Tag ELISA 0.1μg Human Activin A, No Tag Per Well

3.0 Wean Abs. 2.0-0.5
3.0
3.0
2.5
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0
3.0

Log Human Activin RIIB, His Tag Conc.(μg/ml)

-4

-3


Immobilized Human Activin A, No Tag at $1\mu g/ml$ (100 $\mu l/well$) on the plate. Dose response curve for Human Activin RIIB, His Tag with the EC50 of 21.4ng/ml determined by ELISA.

KAGTUS

Assay Data

ELISA Data


Human Activin A, No Tag ELISA 0.5μg Human Activin A, No Tag Per Well



Immobilized Human Activin A, No Tag at $5\mu g/ml$ (100 $\mu l/well$) on the plate. Dose response curve for Human Activin RIIA, His Tag with the EC50 of 0.23 $\mu g/ml$ determined by ELISA.

Cell Based Assay

Recombinant Human Activin A Bioactivity



Measured by its ability to inhibit proliferation of MPC-11 cells. The ED50 for this effect is < 7 ng/mL (QC Test).