

Mouse ADAM8 Protein



Cat. No. ADM-MM108


Description	
Source	Recombinant Mouse ADAM8 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Pro17-Ser658.
Accession	Q05910
Molecular Weight	The protein has a predicted MW of 71.7 kDa. The protein migrates to 55-70 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Supplied as 0.22µm filtered solution in 20mM Tris, 5mM CaCl2, 75mM NaCl, 50% Glycerol (pH 7.5).
Storage	Valid for 12 months from date of receipt when stored at -80°C.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
A Disintegrin And Metalloproteinase (ADAM) proteases constitute a family of multifunctional, membrane-bound proteins with traditional sheddase functions. Their protumorigenic potential has been attributed to both, essential (ADAM10 and ADAM17) and 'dispensable' ADAM proteases (ADAM8, 9, 12, 15, and 19). Of specific interest in this review is the ADAM proteinase ADAM8 that has been identified as a significant player in aggressive malignancies including breast, pancreatic, and brain cancer.	

Assay Data

Bis-Tris PAGE



Mouse ADAM8 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Mouse ADAM8 Protein

Cat. No. ADM-MM108



Assay Data



The purity of Mouse ADAM8 is greater than 95% as determined by SEC-HPLC.

Bioactivity Data

Measured by its ability to cleave a fluorogenic peptide substrate Mca-PLAQAV-Dpa-RSSSR-NH₂. The specific activity is > 0.4 pmol/min/μg.