

Human ANGPTL3/Angiopoietin-like 3 Protein

Cat. No. ANG-HM403

Description

Source	Recombinant Human ANGPTL3/Angiopoietin-like 3 Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Ser17-Lys219.
Accession	Q9Y5C1
Molecular Weight	The protein has a predicted MW of 26.6 kDa. Due to glycosylation, the protein migrates to 30-38 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

ANGPTL3 is a secreted glycoprotein that is structurally related to the angiopoietins. Mature human ANGPTL3 contains an N-terminal coiled coil domain and a C-terminal fibrinogen-like domain. ANGPTL3 is expressed in the liver from early in development through adulthood. Acts in part as a hepatokine that is involved in regulation of lipid and glucose metabolism. Proposed to play a role in the trafficking of energy substrates to either storage or oxidative tissues in response to food intake.

Assay Data

Bis-Tris PAGE




Human ANGPTL3 on Bis-Tris PAGE under reduced conditions. The purity is greater than 95%.

ELISA Data

Human ANGPTL3, His Tag ELISA


0.1µg Human ANGPTL3, His Tag Per Well



Immobilized Human ANGPTL3, His Tag at 1 µg/ml (100µl/Well). Dose response curve for Anti-ANGPTL3 Antibody, hFc Tag with the EC50 of 42.3ng/ml determined by ELISA (QC Test).

Assay Data

SPR Data



Human LILRB2 captured on Protein A chip, can bind Human ANGPTL3, His Tag with an affinity constant of 17.6nM as determined in a SPR assay (Biacore T200).