

# Mouse CD13/ANPEP Protein, Ultra Low Endotoxin

Cat. No. ANP-MM101-UL

## Description

Source	Recombinant Mouse CD13/ANPEP Protein is expressed from HEK293 with His tag at the C-Terminus.
	It contains Lys69-Ser966.
Accession	P97449
Molecular Weight	The protein has a predicted MW of 103.9 kDa. Due to glycosylation, the protein migrates to 115-130 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	> 95% as determined by HPLC

## Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

CD13/aminopeptidase N is a widely expressed ectoenzyme with multiple functions. As an enzyme, CD13 regulates activities of numerous cytokines by cleaving their N-terminals and is involved in Ag processing by trimming the peptides bound to MHC class II. Independent of its enzymatic activity, cell membrane CD13 functions by cross-linking-induced signal transduction, regulation of receptor recycling, enhancement of FcγR-mediated phagocytosis, and acting as a receptor for cytokines.


## Assay Data

### Bis-Tris PAGE



Mouse CD13 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### SEC-HPLC



The purity of Mouse CD13 is greater than 95% as determined by SEC-HPLC.

## Mouse CD13/ANPEP Protein, Ultra Low Endotoxin

Cat. No. ANP-MM101-UL



### Assay Data

#### Bioactivity Data

Measured by its ability to cleave the fluorogenic peptide substrate, Ala-7-amido-4-methylcoumarin (Ala-AMC). The specific activity is >1,000 pmol/min/μg.