

Mouse ANTXR2 Protein

Cat. No. ANT-MM1R2

Description

Source	Recombinant Mouse ANTXR2 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Gln32-Gly318.
Accession	NP_598499.1
Molecular Weight	The protein has a predicted MW of 32.00 kDa. Due to glycosylation, the protein migrates to 33-38 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C . Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

The Capillary Morphogenesis Gene 2 (CMG2) gene encodes an Anthrax toxin receptor (ANTXR2), ANTXR2/CMG2 was originally identified as a result of up-regulation during capillary morphogenesis of endothelial cells (ECs) cultured in vitro. ANTXR2/CMG2 functions to promote endothelial proliferation and morphogenesis during sprouting angiogenesis, consistent with the endothelial expression of ANTXR2/CMG2 in several vascular beds.

Assay Data


Bis-Tris PAGE



Mouse ANTXR2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse ANTXR2 is greater than 95% as determined by SEC-HPLC.