

# Human APLN Protein, Ultra Low Endotoxin

Cat. No. APN-HM201-UL

Description	
Source	Recombinant Human APLN Protein is expressed from HEK293 with hFc (IgG1) tag at the C-Terminus. It contains Gly23-Phe77.
Accession	Q9ULZ1
Molecular Weight	The protein has a predicted MW of 32.9 kDa. Due to glycosylation, the protein migrates to 35-40 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS, 200mM Arginine (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 24 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

**Background**

Macrophages play key roles during cardiovascular diseases (CVD) and their related complications. Apelin (APLN) is a key molecule, whose roles during CVD have been documented previously. Therefore, it has been hypothesized that APLN may perform its roles via modulation of macrophages. Additionally, due to the widespread distribution of the CVD, more effective therapeutic strategies need to be developed to overcome the related complications.

## Assay Data

### Bis-Tris PAGE




Human APLN on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### ELISA Data

Assay Data

**Human APLN, hFc Tag ELISA**

0.1µg Human APLN, hFc Tag Per Well



Immobilized Human APLN, hFc Tag at 1µg/ml (100µl/Well) on the plate. Dose response curve for Biotinylated Anti-APLN Antibody, hFc Tag with the EC50 of 42.2ng/ml determined by ELISA.

Log Biotinylated Anti-APLN Antibody, hFc Tag Conc.(µg/ml)