

Mouse ANGPT2/Angiopoietin-2 Protein

Cat. No. APN-MM101

Description

Source	Recombinant Mouse ANGPT2/Angiopoietin-2 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Lys275-Phe496.
Accession	O35608
Molecular Weight	The protein has a predicted MW of 26.29 kDa. Due to glycosylation, the protein migrates to 30-35 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Angiopoietin-2 (Ang-2; also ANGPT2) is a secreted glycoprotein that plays a complex role in angiogenesis and inflammation. Mature Ang-2 is 478 amino acids in length. Ang2 is widely expressed during development, but it is restricted postnatally to highly angiogenic tissues such as the placenta, ovaries, and uterus. It is particularly abundant in vascular endothelial cells (EC) where it is stored in intracellular Weibel Palade bodies..

Assay Data

Bis-Tris PAGE



Mouse ANGPT2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Mouse ANGPT2/Angiopoietin-2 Protein

Cat. No. APN-MM101



Assay Data



The purity of Mouse ANGPT2 is greater than 95% as determined by SEC-HPLC.