

Human APOE3/Apolipoprotein E Protein



Cat. No. APO-HM101


Description	
Source	Recombinant Human APOE3/Apolipoprotein E Protein is expressed from HEK293 with His tag at the N-Terminus. It contains Lys19-His317.
Accession	P02649-1
Molecular Weight	The protein has a predicted MW of 35.3 kDa. Due to glycosylation, the protein migrates to 35-40 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage	
Formulation	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Apolipoprotein E (apoE) is a lipid carrier in both the peripheral and the central nervous systems. Lipid-loaded apoE lipoprotein particles bind to several cell surface receptors to support membrane homeostasis and injury repair in the brain. Considering prevalence and relative risk magnitude, the ε4 allele of the APOE gene is the strongest genetic risk factor for late-onset Alzheimer's disease (AD).	

Assay Data

Bis-Tris PAGE



Human APOE3 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.