

Human APOE4/Apolipoprotein E Protein, Ultra Low Endotoxin

Cat. No. APO-HM202-UL

Description

Source	Recombinant Human APOE4/Apolipoprotein E Protein is expressed from HEK293 with hFc (IgG1) tag at the N-Terminus. It contains Lys19-His317.
Accession	AAB59397.1
Molecular Weight	The protein has a predicted MW of 61.6 kDa. Due to glycosylation, the protein migrates to 63-67 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage

Formulation	Lyophilized from 0.22µm filtered solution in PBS, 2mM DTT, 5mM CHAPS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Apolipoprotein E (apoE) is a lipid carrier in both the peripheral and the central nervous systems. Lipid-loaded apoE lipoprotein particles bind to several cell surface receptors to support membrane homeostasis and injury repair in the brain. Considering prevalence and relative risk magnitude, the ε4 allele of the APOE gene is the strongest genetic risk factor for late-onset Alzheimer's disease (AD).

Assay Data

Bis-Tris PAGE



Human APOE4 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.