
## Mouse APOE/Apolipoprotein E Protein





Description	
Source	Recombinant Mouse APOE/Apolipoprotein E Protein is expressed from HEK293 with hFc tag at the C-Terminus.
	It contains Glu19-Gln311.
Accession	P08226
Molecular Weight	The protein has a predicted MW of 60.7 kDa. Due to glycosylation, the protein migrates to 62-66 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC

#### Formulation and Storage


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in 20mM PB, 150mM NaCl (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt20 to -80°C for 3-6 months in unopened state after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

# **Background**

Apolipoprotein E (apoE) is a lipid carrier in both the peripheral and the central nervous systems. Lipid-loaded apoE lipoprotein particles bind to several cell surface receptors to support membrane homeostasis and injury repair in the brain. Considering prevalence and relative risk magnitude, the ε4 allele of the APOE gene is the strongest genetic risk factor for late-onset Alzheimer's disease (AD).

## **Assay Data**

### **Tris-Bis PAGE**




Mouse APOE on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Cat. No. APO-MM202



# **Assay Data**



The purity of Mouse APOE is greater than 95% as determined by SEC-HPLC.