Biotinylated Human B7-2/CD86 Protein


Cat. No. B72-HM486B

K∧₲℃℧

Description	
Source	Recombinant Biotinylated Human B7-2/CD86 Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus.
	It contains Leu26-Pro247.
Accession	P42081-1
Molecular Weight	The protein has a predicted MW of 28.2 kDa. Due to glycosylation, the protein migrates to 55-70 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μ g by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt80°C for 3-6 months after reconstitution. 2-8°C for 2- 7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Background	
	B7-1 and B7-2 are homologous costimulatory ligands expressed on the surface of antigen presenting cells (APCs). Binding of these molecules to the T cell costimulatory receptors, CD28 and CTLA-4, is essential for the activation and regulation of T cell immunity. B7-1 and B7-2 do not form hetero-oligomers, underscoring the biological relevance of dimeric and monomeric state of B7-1 and B7-2, respectively.

Assay Data


Tris-Bis PAGE



Biotinylated Human B7-2 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

ELISA Data





Immobilized Biotinylated Human B7-2, His Tag at 1μ g/ml (100 μ l/well) on the streptavidin precoated plate (5 μ g/ml). Dose response curve for Human CTLA-4, hFc Tag with the EC50 of 1.9ng/ml determined by ELISA.