

Human BAFF/TNFSF13B/CD257 (68-133) Protein

Cat. No. BAF-HM13D

Description

Source	Recombinant Human BAFF/TNFSF13B/CD257 (68-133) Protein is expressed from HEK293 with His tag at the N-terminus. It contains Gln68-Arg133.
Accession	Q9Y275-1
Molecular Weight	The protein has a predicted MW of 7.77 kDa. Due to glycosylation, the protein migrates to 20-30 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 90% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

B-cell activating factor (BAFF) also known as tumor necrosis factor ligand superfamily member 13B is a protein that in humans is encoded by the TNFSF13B gene. BAFF is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B/TACI, TNFRSF17/BCMA, and TNFRSF13C/BAFF-R.


Assay Data

Bis-Tris PAGE



Human BAFF (68-133) on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Human BAFF (68-133) is greater than 90% as determined by SEC-HPLC.