FITC-Compatible Human BCMA/TNFRSF17 Protein






Description	
Source	Recombinant FITC-Compatible Human BCMA/TNFRSF17 Protein is expressed from HEK293 with His tag at the C-Terminus.
	It contains Met1-Ala54.
Accession	Q02223-1
Molecular Weight	The protein has a predicted MW of 35.9 kDa. Due to glycosylation, the protein migrates to 37-48 kDa based on Bis-Tris PAGE result.
Wavelength	Excitation Wavelength: 490 nm
	Emission Wavelength: 520 nm
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	> 95% as determined by HPLC
Formulation and Storage	
Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Background	
	B-cell maturation antigen (BCMA or BCM), also known as tumor necrosis factor receptor superfamily member 17 (TNFRSF17), is a protein that in humans is encoded by the TNFRSF17 gene.TNFRSF17 is a cell surface receptor of the TNF receptor superfamily which recognizes B-cell activating factor (BAFF).

Assay Data

Bis-Tris PAGE




FITC-Compatible Human BCMA on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Cat. No. BCM-HM117C





Assay Data




The purity of FITC-Compatible Human BCMA is greater than 95% as determined by SEC-HPLC.

FACS Data







Use 100 μ I FITC-Compatible Human BCMA (1 μ g/ml) to detect the positive rate of 1×10⁶ anti-BCMA CAR cells. Non-transfected 293T cells and FITC-Labeled irrelevant protein (100 μ I, 1 μ g/ml) were served as negative control.