

# Mouse B7-H5/Gi24/VISTA Protein

Cat. No. BH7-MM275



## Description

<b>Source</b>	Recombinant Mouse B7-H5/Gi24/VISTA Protein is expressed from HEK293 with hFc tag at the C-Terminus. It contains Phe33-Ala194.
<b>Accession</b>	NP_083008.1
<b>Molecular Weight</b>	The protein has a predicted MW of 44.77 kDa. Due to glycosylation, the protein migrates to 60-68 kDa based on Bis-Tris PAGE result.
<b>Endotoxin</b>	Less than 1EU per µg by the LAL method.
<b>Purity</b>	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

## Formulation and Storage


<b>Formulation</b>	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
<b>Storage</b>	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

B7-H5, also known as VISTA, B7H5, Dies1, SISP1 and C10orf54, is a 55-65 kDa member of the Ig superfamily. It is a transmembrane molecule expressed in bone, on embryonic stem cells (ESCs), and on tumor cell surfaces.


## Assay Data

### Bis-Tris PAGE



Mouse B7-H5 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

### SEC-HPLC



The purity of Mouse B7-H5 is greater than 95% as determined by SEC-HPLC.