

Mouse CA9/Carbonic Anhydrase IX Protein, Ultra Low Endotoxin

Cat. No. CA9-MM101-UL

Description

Source	Recombinant Mouse CA9/Carbonic Anhydrase IX Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Gln32-Asp390.
Accession	NP_647466.2
Molecular Weight	The protein has a predicted MW of 43.12 kDa. Due to glycosylation, the protein migrates to 50-53 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CA9 is a member of the carbonic anhydrases' family, that is often expressed in cancer cells under hypoxic condition. CA9 expression potentially contributes to the regulation of cancer cell differentiation and mediates tumour-associated genes and signalling pathways, including apoptosis, hypoxia, G2M checkpoint, PI3K/AKR/mTOR signalling and TGF-beta signalling pathways.

Assay Data


Bis-Tris PAGE



Mouse CA9 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse CA9 is greater than 95% as determined by SEC-HPLC.