

Cynomolgus Calprotectin (S100A8&S100A9) Protein

Cat. No. CAL-CE101



Description

Source	Recombinant Cynomolgus Calprotectin (S100A8&S100A9) Protein is expressed from E.coli with His tag at the C-terminus. It contains Met1-Glu93 (S100A8) and Met1-Arg114 (S100A9).
Accession	A0A2K5X560(S100A8)&XP_005541846.1(S100A9)
Molecular Weight	The protein has a predicted MW of 11.91 kDa (S100A8) and 13.47 kDa (S100A9) same as Bis-Tris PAGE result.
Endotoxin	Less than 1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 90% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Calprotectin, a member of the widespread calcium-binding S-100 protein family, is present in remarkably high concentration in the cytoplasm of human neutrophils. The calprotectin molecule is composed of light (MRP8) and heavy (MRP14) subunits. Although not secreted from intact neutrophils, calprotectin release from dead and dying neutrophils creates high concentrations of the protein in inflammatory or abscess fluids and in the intestinal tract lumen of patients with inflammatory bowel disease.


Assay Data

Bis-Tris PAGE



Cynomolgus Calprotectin (S100A8&S100A9) on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Cynomolgus Calprotectin (S100A8&S100A9) is greater than 90% as determined by SEC-HPLC.