

Human CALCA/CGRP Protein, Ultra Low Endotoxin

Cat. No. CAL-HM201-UL

Description

Source	Recombinant Human CALCA/CGRP Protein is expressed from HEK293 with hFc (IgG1) tag at the C-Terminus. It contains Ala26-Asn141.
Accession	NP_001732.1
Molecular Weight	The protein has a predicted MW of 39.4 kDa. Due to glycosylation, the protein migrates to 45-48 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Specific fragments of methylation changes in the target gene (Calca) were identified by IGV analysis. CGRP was applied to compare the effects on ASCs-T2DM morphology via phalloidin staining, proliferation through CCK-8 assay, and osteogenic differentiation with osteogenic staining, qPCR, and repair of calvarial defect. Furthermore, 5-azacytidine (5-az) was used to intervene ASCs-T2DM to verify the relationship between the methylation level of the target fragment and expression of Calca.

Assay Data

Bis-Tris PAGE



Human CALCA on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.