

Cynomolgus PTK7/CCK4 Protein, Ultra Low Endotoxin



Cat. No. CCK-CM104-UL

Description	
Source	Recombinant Cynomolgus PTK7/CCK4 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Ala39-Thr712.
Accession	A0A2K5VNI5
Molecular Weight	The protein has a predicted MW of 75.73 kDa. Due to glycosylation, the protein migrates to 90-110 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Protein Tyrosine Kinase 7 (PTK7) is as a critical regulator of canonical and non-canonical Wnt-signaling during embryonic development and cancer cell formation. Disrupting PTK7 activity perturbs vertebrate nervous system development, and also promotes human cancer formation. Observations in different model systems suggest a complex cross-talk between PTK7 protein and Wnt signaling.	

Assay Data


Bis-Tris PAGE



Cynomolgus PTK7 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Cynomolgus PTK7 is greater than 95% as determined by SEC-HPLC.