

Non-biotinylated Human CD160 Protein, Ultra Low Endotoxin

Cat. No. CD1-HM460-UL

Description

Source	Recombinant Human CD160 Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Gly25-Leu158.
Accession	O95971-1
Molecular Weight	The protein has a predicted MW of 17.7 kDa. Due to glycosylation, the protein migrates to 27-33 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD160 (also Natural killer cell receptor BY55) is a 27-30 kDa member of the Ig superfamily. In human, it is expressed principally on nonmyeloid hematopoietic cells. CD160 antigen is a receptor on immune cells capable to deliver stimulatory or inhibitory signals that regulate cell activation and differentiation. Exists as a GPI-anchored and as a transmembrane form, each likely initiating distinct signaling pathways via phosphoinositol 3-kinase in activated NK cells and via LCK and CD247/CD3 zeta chain in activated T cells.

Assay Data


Bis-Tris PAGE



Human CD160 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Assay Data




The purity of Human CD160 is greater than 95% as determined by SEC-HPLC.

BLI Data



Loaded Human HVEM, hFc Tag on ProA-Biosensor can bind Human CD160, His Tag with an affinity constant of $0.95 \mu\text{M}$ as determined in BLI assay .

SPR Data



Human HVEM, hFc Tag captured on CM5 Chip via Protein A can bind Human CD160, His Tag with an affinity constant of $1.44 \mu\text{M}$ as determined in SPR assay (Biacore T200).