

Mouse CD180 Protein

Cat. No. CD1-MM180

Description

Source	Recombinant Mouse CD180 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Thr21-Ser626.
Accession	Q62192
Molecular Weight	The protein has a predicted MW of 69.1 kDa. Due to glycosylation, the protein migrates to 80-110 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C . Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD180, a related member of the Toll-like receptor family, is lost or underexpressed at the plasma membrane in circulating cells of various B-cell lymphomas except marginal zone lymphomas (MZL). CD180 can be considered as a single positive robust marker of MZL and should be therefore included in flow cytometry panels for the diagnosis of mature B-cell neoplasms.

Assay Data


Bis-Tris PAGE



Mouse CD180 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse CD180 is greater than 95% as determined by SEC-HPLC.