

Cynomolgus CD24 Protein

Cat. No. CD2-CM224



Description

Source	Recombinant Cynomolgus CD24 Protein is expressed from Expi293 with hFc tag at the C-terminal. It contains Thr27-Lys54.
Accession	I7GKK1-1
Molecular Weight	The protein has a predicted MW of 30.2 kDa. Due to glycosylation, the protein migrates to 45-60 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22µm filtered solution in PBS (pH 7.4). Please dilute to the desired concentration according to the concentration of the solution shown on the product label.
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please do not repeated freeze-thaw cycles.

Background

CD24 is a sialoglycoprotein expressed at the surface of most B lymphocytes and differentiating neuroblasts. It is also expressed on neutrophils and neutrophil precursors from the myelocyte stage onwards. The potential for targeting CD24 in cancer therapy seems promising, as CD24 is overexpressed in many human cancers.


Assay Data

Tris-Bis PAGE



Cynomolgus CD24 on Tris-Bis PAGE under reduced conditions. The purity is greater than 95%.

SEC-HPLC



The purity of Cynomolgus CD24 is greater than 95% as determined by SEC-HPLC.