Human CD20 Protein-VLP

Cat. No. CD2-HM122



Description	
Source	Recombinant Human CD20 Protein-VLP is expressed from HEK293.
	It contains Met1-Pro297.
Accession	P11836-1
Molecular Weight	The target protein has a predicted MW of 34.3 kDa.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by HPLC

Formulation and Storage

Formulation	Supplied as 0.22µm filtered solution in PBS, 300mM L-Arginine (pH 7.4). Notice: If you need it for immunization,
	Do Not use any adjuvant


Storage Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

B-lymphocyte antigen CD20 or CD20 is an activated-glycosylated phosphoprotein expressed on the surface of all B-cells beginning at the pro-B phase (CD45R, CD117) and progressively increasing in concentration until maturity.CD20 is the target of the monoclonal antibodies rituximab, ocrelizumab, obinutuzumab, ofatumumab, ibritumomab tiuxetan, tositumomab, and ublituximab, which are all active agents in the treatment of all B cell lymphomas, leukemias, and B cell-mediated autoimmune diseases.

Assay Data


SEC-HPLC



The purity of Human CD20 VLP is greater than 95% as determined by SEC-HPLC.

ELISA Data

Human CD20 VLP ELISA 0.5µg Human CD20 VLP Per Well




Immobilized Human CD20 VLP at 5µg/ml (100µl/Well) on plate on the plate. Dose response curve for Rituximab, hFc Tag with the EC50 of 11.4ng/ml determined by ELISA (QC Test).

KAGTUS

Assay Data

ELISA Data


Human CD20 VLP ELISA 0.5µg Human CD20 VLP Per Well



Immobilized Human CD20 VLP at $5\mu g/ml$ (100 $\mu l/Well$) on plate on the plate. Dose response curve for Ofatumumab, hFc Tag with the EC50 of 10.5ng/ml determined by ELISA (QC Test).

ELISA Data

Human CD20 VLP ELISA 0.5µg Human CD20 VLP Per Well



Immobilized Human CD20 VLP at 5µg/ml (100µl/well) on the plate. Dose response curve for Obinutuzumab, hFc Tag with the EC50 of 1.7ng/ml determined by ELISA.