## Human CD24 Protein, Ultra Low Endotoxin

Cat. No. CD2-HM324-UL



| Description         |                                                                                                                                    |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Source              | Recombinant Human CD24 Protein is expressed from HEK293 with mFc (IgG1) tag at the C-Terminus.                                     |
|                     | It contains Ser27-Gly59.                                                                                                           |
| Accession           | P25063-1                                                                                                                           |
| Molecular<br>Weight | The protein has a predicted MW of 29.5 kDa. Due to glycosylation, the protein migrates to 48-60 kDa based on Bis-Tris PAGE result. |
| Endotoxin           | Less than 0.001 EU per μg by the LAL method.                                                                                       |
| Purity              | > 95% as determined by Bis-Tris PAGE                                                                                               |
|                     | > 95% as determined by HPLC                                                                                                        |
| Formulation and     | l Storage                                                                                                                          |


| Formulation    | Lyophilized from 0.22μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.                                                                                      |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reconstitution | Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.                                                                                         |
| Storage        | -20 to -80°C for 12 months as supplied from date of receipt80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles. |
|                |                                                                                                                                                                                                                     |

### Background

CD24 is a sialoglycoprotein expressed at the surface of most B lymphocytes and differentiating neuroblasts. It is also expressed on neutrophils and neutrophil precursors from the myelocyte stage onwards. The potential for targeting CD24 in cancer therapy seems promising, as CD24 is overexpressed in many human cancers.

### **Assay Data**

### **Bis-Tris PAGE**




Human CD24 on Bis-Tris PAGE under reduced conditions. The purity is greater than 95%.

**SEC-HPLC** 

# KAGTUS


### **Assay Data**



The purity of Human CD24 is greater than 95% as determined by SEC-HPLC.

### **ELISA Data**

# Human CD24, mFc Tag ELISA 0.2µg Human CD24, mFc Tag Per Well



Immobilized Human CD24, mFc Tag at 2µg/ml (100µl/well) on the plate. Dose response curve for Anti-CD24 Antibody, hFc Tag with the EC50 of 54.3ng/ml determined by ELISA (QC Test).