

# Human CD36/SR-B3 Protein

Cat. No. CD3-HM236

## Description

<b>Source</b>	Recombinant Human CD36/SR-B3 Protein is expressed from HEK293 with hFc (IgG1) tag at the C-terminus. It contains Gly30-Asn439.
<b>Accession</b>	P16671-1
<b>Molecular Weight</b>	The protein has a predicted MW of 72.56 kDa. Due to glycosylation, the protein migrates to 90-120 kDa based on Bis-Tris PAGE result.
<b>Endotoxin</b>	Less than 0.1 EU per µg by the LAL method.
<b>Purity</b>	> 90% as determined by Bis-Tris PAGE > 95% as determined by HPLC

## Formulation and Storage


<b>Formulation</b>	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
<b>Storage</b>	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

CD36, alternatively known as platelet membrane glycoprotein IV (GPIV), GPIIb, thrombospondin receptor, collagen receptor, fatty acid translocase (FAT), and scavenger receptor class B, member 3 (SR-B3), is an integral membrane glycoprotein that has multiple physiological functions. CD36 is a multifunctional glycoprotein that acts as receptor for a broad range of ligands. Ligands can be of proteinaceous nature like thrombospondin, fibronectin, collagen or amyloid-beta as well as of lipidic nature such as oxidized low-density lipoprotein (oxLDL), anionic phospholipids, long-chain fatty acids and bacterial diacylated lipopeptides.

## Assay Data


### Bis-Tris PAGE



Human CD36 on Bis-Tris PAGE under reduced condition. The purity is greater than 90%.

### SEC-HPLC

Assay Data



The purity of Human CD36 is greater than 95% as determined by SEC-HPLC.