Mouse CD31/PECAM-1 Protein

Cat. No. CD3-MM131



Description	
Source	Recombinant Mouse CD31/PECAM-1 Protein is expressed from HEK293 with His tag at the C-Terminus.
	It contains Glu18-Lys590.
Accession	Q08481-1
Molecular Weight	The protein has a predicted MW of 65.6 kDa. Due to glycosylation, the protein migrates to 90-100 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt20 to -80°C for 3-6 months in unopened state after reconstitution.2-8°C for 2-7 days after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Platelet endothelial cell adhesion molecule 1 (PECAM-1) is an adhesion and signaling receptor that is expressed on endothelial and hematopoietic cells and plays important roles in angiogenesis, vascular permeability, and regulation of cellular responsiveness.

Assay Data

Tris-Bis PAGE




Mouse CD31/PECAM-1 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Cat. No. CD3-MM131



Assay Data



The purity of Mouse CD31/PECAM-1 is greater than 95% as determined by SEC-HPLC.