

Human CD46 Protein

Cat. No. CD4-HM146

Description

Source	Recombinant Human CD46 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Cys35-Asp343.
Accession	P15529-1
Molecular Weight	The protein has a predicted MW of 35.4 kDa. Due to glycosylation, the protein migrates to 60-70 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage

Formulation	Supplied as 0.22µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD46 was discovered in 1986 during a search for novel C3b-binding proteins. CD46 is expressed ubiquitously and functions as a co-factor in the factor I-mediated proteolytic cleavage of C3b and C4b. Its vital role in preventing complement deposition on host tissue is underpinned by the fact that deficiency of CD46 is a predisposing factor for numerous disease conditions arising from complement-mediated 'self-attack'.

Assay Data


Bis-Tris PAGE



Human CD46 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Human CD46 is greater than 95% as determined by SEC-HPLC.

ELISA Data

Human CD46, His Tag ELISA

0.05µg Human CD46, His Tag Per Well



Immobilized Human CD46, His Tag at 0.5µg/ml (100µl/Well) on the plate. Dose response curve for Anti-CD46 Antibody, hFc Tag with the EC50 of 4.2ng/ml determined by ELISA.