

Mouse CD164 Protein

Cat. No. CD6-MM204

Description

Source	Recombinant Mouse CD164 Protein is expressed from HEK293 with hFc tag at the C-Terminus. It contains Gln24-Asp162.
Accession	Q9R0L9
Molecular Weight	The protein has a predicted MW of 41.4 kDa. Due to glycosylation, the protein migrates to 70-115 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -20 to -80°C for 3-6 months in unopened state after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD164 was found to play a role in many malignant diseases. CD164 was associated with clinical and pathological features of patients. High level of CD164 was related to the distant metastasis and vascular invasion of bladder cancer patients. CD164 was associated with the poor clinical outcomes of BC patients. Silencing of CD164 could inhibit the progression of tumors in vivo and in vitro, which may become an effective target in the treatment of bladder cancer.

Assay Data

Tris-Bis PAGE




Mouse CD164 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Mouse CD164 Protein

Cat. No. CD6-MM204

Assay Data



The purity of Mouse CD164 is greater than 95% as determined by SEC-HPLC.