

Human CD79A Protein, Ultra Low Endotoxin

Cat. No. CD7-HM19A-UL



Description

Source	Recombinant Human CD79A Protein is expressed from HEK293 with His tag at the C-terminus.
	It contains Leu33-Arg143.
Accession	P11912-1
Molecular Weight	The protein has a predicted MW of 14.17 kDa. Due to glycosylation, the protein migrates to 40-60 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	> 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD79 is a heterodimeric molecule comprising two polypeptide chains, B29 (CD79b) and mb-1 (CD79a). It is physically linked in the surface of B cells to membrane immunoglobulin, forming the B cell antigen receptor complex. Expression of the mb-1 (CD79a) chain has been studied in leukaemias and shown to be present in most B lineage acute lymphoblastic leukaemias (ALL).

Assay Data

Bis-Tris PAGE



Human CD79A on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Human CD79A Protein, Ultra Low Endotoxin

Cat. No. CD7-HM19A-UL



Assay Data



The purity of Human CD79A is greater than 95% as determined by SEC-HPLC.