Mouse CD27 Ligand/CD70 Protein


Cat. No. CD7-MM270

ϗͶϲϿ·Ͷϩ

Description	
Source	Recombinant Mouse CD27 Ligand/CD70 Protein is expressed from HEK293 with hFc tag at the N-Terminus.
	It contains Ser45-Pro195.
Accession	O55237
Molecular Weight	The protein has a predicted MW of 43.41 kDa. Due to glycosylation, the protein migrates to 55-60 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	> 95% as determined by HPLC
Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 6.2). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 24 months as supplied from date of receipt80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.
Background	
	CD70, also named CD27 ligand (CD27L), is a type II transmembrane glycoprotein belonging to the TNF superfamily (TNFSF) and has been designated TNFSF7. CD70 is a cytokine that binds to CD27. Plays a role in T-cell activation. Induces the proliferation of costimulated T-cells and enhances the generation of cytolytic T-cells.


Assay Data

Bis-Tris PAGE



Mouse CD27 Ligand on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Mouse CD27 Ligand is greater than 95% as determined by SEC-HPLC.

Mouse CD27 Ligand/CD70 Protein


Cat. No. CD7-MM270

Assay Data

ELISA Data

ϗͺϒͺϿ·ϽϨ





Immobilized Mouse CD27, mFc Tag at 2μ g/ml (100 μ l/Well) on the plate. Dose response curve for Mouse CD27 Ligand, hFc Tag with the EC50 of 34.0ng/ml determined by ELISA.