

Biotinylated Human CD94 Protein

Cat. No. CD9-HM404B

Description

Source	Recombinant Biotinylated Human CD94 Protein is expressed from HEK293 with His tag and Avi tag at the N-Terminus. It contains Ser34-Ile179.
Accession	Q13241-1
Molecular Weight	The protein has a predicted MW of 19.8 kDa. Due to glycosylation, the protein migrates to 35-45 kDa based on Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -20 to -80°C for 3-6 months in unopened state after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD94 is an approximately 25 kDa type 2 transmembrane protein that plays an important role in regulating natural killer (NK) cell activation. CD94 plays a role as a receptor for the recognition of MHC class I HLA-E molecules by NK cells and some cytotoxic T-cells.

Assay Data

Tris-Bis PAGE



Biotinylated Human CD94 on Tris-Bis PAGE under reduced conditions. The purity is greater than 95%.


SEC-HPLC

Biotinylated Human CD94 Protein

Cat. No. CD9-HM404B



Assay Data



The purity of Biotinylated Human CD94 is greater than 95% as determined by SEC-HPLC.