

Human Fc gamma RIIB/C (CD32b/c) Protein, Ultra Low Endotoxin



Cat. No. CDB-HM401-UL


Description	
Source	Recombinant Human Fc gamma RIIB/C (CD32b/c) Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus. It contains Ala46-Pro217.
Accession	P31994-1
Molecular Weight	The protein has a predicted MW of 22.5 kDa. Due to glycosylation, the protein migrates to 32-37 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt.-80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
The Fc gamma Rs have been divided into three classes based on close relationships in their extracellular domains; these groups are designated Fc gamma RI (also known as CD64), Fc gamma RII (CD32), and Fc gamma RIII (CD16). Each group may be encoded by multiple genes and exist in different isoforms depending on species and cell type.	

Assay Data


Bis-Tris PAGE



Human Fc gamma RIIB/C on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Assay Data




The purity of Human Fc gamma RIIB/C is greater than 95% as determined by SEC-HPLC.

SPR Data



Human Fc gamma RIIB/C, His-Avi Tag captured on CM5 Chip via Anti-his antibody can bind Rituximab with an affinity constant of 16.72 μM as determined in SPR assay (Biacore T200) (QC Test).



Rituximab captured on CM5 Chip via Protein A can bind Human Fc gamma RIIB/C, His Tag with an affinity constant of 3.46 μM as determined in SPR assay (Biacore T200).