

Human CDH17/Cadherin 17 Domain 7 Protein

Cat. No. CDH-HM37D

Description

Source	Recombinant Human CDH17/Cadherin 17 Domain 7 Protein is expressed from HEK293 with mFc (IgG2a) tag at the C-terminus. It contains Ala668-Met787.
Accession	Q12864
Molecular Weight	The protein has a predicted MW of 39.54 kDa. Due to glycosylation, the protein migrates to 48-58 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.1EU per µg by the LAL method.
Purity	>95% as determined by Bis-Tris PAGE >95% as determined by HPLC

Formulation and Storage

Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Liver-intestine cadherin (CDH17) has been known to function as a tumor stimulator and diagnostic marker for almost two decades. In vivo studies showed CDH17 knockout resulted in apoptotic PC tumor death through activating caspase-3 activity. Taken together, CDH17 functions as an oncogenic molecule critical to PC growth by regulating tumor apoptosis signaling pathways and CDH17 could be targeted to develop an anti-PC therapeutic approach.


Assay Data

Bis-Tris PAGE



Human CDH17 Domain 7 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Human CDH17 Domain 7 is greater than 95% as determined by SEC-HPLC.