

Biotinylated Human CD40 Ligand/TNFSF5 Trimer Protein (Primary Amine Labeling)



Cat. No. CDL-HM140B

Description

Source	Recombinant Biotinylated Human CD40 Ligand/TNFSF5 Trimer Protein (Primary Amine Labeling) is expressed from HEK293 with His tag and Flag tag at the N-Terminus. It contains Met113-Leu261.
Accession	P29965
Molecular Weight	The protein has a predicted MW of 52.4 kDa, Due to glycosylation, the protein migrates to 55-60 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD40 ligand or CD40L, also called CD154, is a protein that is primarily expressed on activated T cells and is a member of the TNF superfamily of molecules. It binds to CD40 (protein) on antigen-presenting cells (APC), which leads to many effects depending on the target cell type. In total CD40L has three binding partners: CD40, $\alpha 5\beta 1$ integrin and $\alpha 1\text{Ib}\beta 3$. CD154 acts as a costimulatory molecule and is particularly important on a subset of T cells called T follicular helper cells.


Assay Data

Bis-Tris PAGE



Biotinylated Human CD40 Ligand Trimer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC




The purity of Biotinylated Human CD40 Ligand Trimer is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Biotinylated Human CD40 Ligand Trimer, His Tag ELISA

0.5µg Human CD40, hFc Tag Per Well



Immobilized Human CD40, hFc Tag at 5µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Human CD40 Ligand Trimer, His Tag with the EC50 of 47.6ng/ml determined by ELISA (QC Test).

Log Biotinylated Human CD40 Ligand Trimer, His Tag Conc.(µg/ml)