

Mouse CHODL Protein

Cat. No. CDL-MM201

Description

Source	Recombinant Mouse CHODL Protein is expressed from HEK293 with hFc tag at the C-Terminus. It contains Arg22-Asn216.
Accession	Q9CXM0-1
Molecular Weight	The protein has a predicted MW of 48.7 kDa. Due to glycosylation, the protein migrates to 65-68 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

chondrolectin (CHODL) is commonly overexpressed in the majority of lung cancers, indicating a possible correlation between CHODL and metastasis of lung cancer cells. the expression of CHODL is significantly decreased in HCC clinical samples and in HCC cell lines. Overexpression of CHODL in SMMC7721 cells with a lentiviral vector increased SMMC7721 cell migration and invasion.

Assay Data


Bis-Tris PAGE



Mouse CHODL on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse CHODL is greater than 95% as determined by SEC-HPLC.