

Biotinylated Mouse MDL-1/CLEC5A Protein

Cat. No. CLE-MM45AB

Description

Source	Recombinant Biotinylated Mouse MDL-1/CLEC5A Protein is expressed from HEK293 with His tag and Avi tag at the N-terminus. It contains Tyr26-Lys190.
Accession	Q9R007-3
Molecular Weight	The protein has a predicted MW of 21.76 kDa. Due to glycosylation, the protein migrates to 35-50 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μ g by the LAL method.
Purity	>95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Supplied as 0.22 μ m filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

In *C. elegans*, increased lifespan in daf-2 insulin/IGF-1 receptor mutants is accompanied by up-regulation of the MDL-1 Mad basic helix-loop-helix leucine zipper transcription factor. MDL-1, like its mammalian orthologs, is an inhibitor of cell proliferation and growth that slows progression of an age-related pathology in *C. elegans* (uterine tumors). In addition, intestine-limited rescue of mdl-1 increased lifespan but not to wild type levels. Thus, mdl-1 likely acts both in the intestine and the germline to influence age-related mortality.

Assay Data

Bis-Tris PAGE



Biotinylated Mouse MDL-1 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.