

Cynomolgus M-CSF/CSF-1 Protein

Cat. No. CSF-CM101

Description

Source	Recombinant Cynomolgus M-CSF/CSF-1 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Glu33-Leu255.
Accession	XP_045221641.1
Molecular Weight	The protein has a predicted MW of 26.26 kDa. Due to glycosylation, the protein migrates to 45-55 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Supplied as 0.22 μm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C . Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CSF-1, also known as M-CSF, is a four α -helical-bundle cytokine that is the primary regulator of macrophage survival, proliferation and differentiation. CSF-1 is also essential for the survival and proliferation of osteoclast progenitors. CSF-1 is cytokine that plays an essential role in the regulation of survival, proliferation and differentiation of hematopoietic precursor cells, especially mononuclear phagocytes, such as macrophages and monocytes. Promotes the release of proinflammatory chemokines, and thereby plays an important role in innate immunity and in inflammatory processes. This mutation leads to a deletion of exons 2-7 of the EGFR gene and renders the mutant receptor incapable of binding any known ligand.

Assay Data


Bis-Tris PAGE



Cynomolgus M-CSF on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


ELISA Data

Assay Data

Cynomolgus M-CSF, His Tag ELISA0.1 μ g Cynomolgus M-CSF, His Tag Per Well

Immobilized Cynomolgus M-CSF, His Tag at 1 μ g/ml (100 μ l/well) on the plate. Dose response curve for Human M-CSF R, hFc Tag with the EC50 of 6.5ng/ml determined by ELISA (QC Test).

SPR Data



Cynomolgus M-CSF R, His Tag immobilized on CM5 Chip can bind Cynomolgus M-CSF, His Tag with an affinity constant of 0.889 μ M as determined in SPR assay (Biacore T200).