

Mouse CTLA-4/CD152 Protein

Cat. No. CTL-MM141



Description

Source	Recombinant Mouse CTLA-4/CD152 Protein is expressed from HEK293 with His tag at the C-terminus. It contains Glu36-Asp161.
Accession	P09793
Molecular Weight	The protein has a predicted MW of 15.47 kDa. Due to glycosylation, the protein migrates to 28-38 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.1 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Supplied as 0.22 µm filtered solution in PBS (pH 7.4).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CTLA-4 (cytotoxic T-lymphocyte-associated protein 4), also known as CD152, is a protein receptor that, functioning as an immune checkpoint, downregulates immune responses. CTLA4 is constitutively expressed in regulatory T cells but only upregulated in conventional T cells after activation – a phenomenon which is particularly notable in cancers. It acts as an "off" switch when bound to CD80 or CD86 on the surface of antigen-presenting cells.


Assay Data

Bis-Tris PAGE



Mouse CTLA-4 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC




The purity of Mouse CTLA-4 is greater than 95% as determined by SEC-HPLC.

Assay Data

ELISA Data

Mouse CTLA-4, His Tag ELISA

0.1µg Mouse CTLA-4, His Tag Per Well



Immobilized Mouse CTLA-4, His Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Human B7-1, hFc Tag (Cat. B71-HM280) with the EC50 of 1.9ng/ml determined by ELISA.