

# Human CX3CL1/Fractalkine Protein

Cat. No. CX3-HM101

## Description

<b>Source</b>	Recombinant Human CX3CL1/Fractalkine Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Gln25-Gln341.
<b>Accession</b>	AAH01163
<b>Molecular Weight</b>	The protein has a predicted MW of 34.7 kDa. Due to glycosylation, the protein migrates to 90-110 kDa based on Bis-Tris PAGE result.
<b>Endotoxin</b>	Less than 1EU per $\mu$ g by the LAL method.
<b>Purity</b>	> 95% as determined by Bis-Tris PAGE

## Formulation and Storage


<b>Formulation</b>	Supplied as 0.22 $\mu$ m filtered solution in PBS (pH 7.4).
<b>Storage</b>	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

Fractalkine/CX3C chemokine ligand 1 (CX3CL1) is a chemokine involved in the anticancer function of lymphocytes-mainly NK cells, T cells and dendritic cells. Its increased levels in tumors improve the prognosis for cancer patients, although it is also associated with a poorer prognosis in some types of cancers, such as pancreatic ductal adenocarcinoma.

## Assay Data

### Bis-Tris PAGE



Human CX3CL1 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.