

Cynomolgus DDT Protein

Cat. No. DDT-CE101

Description

Source	Recombinant Cynomolgus DDT Protein is expressed from E.coli with His tag at the N-Terminus. It contains Pro2-Leu118.
Accession	XP_005595435.2
Molecular Weight	The protein has a predicted MW of 13.37 kDa same as Tris-Bis PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Tris-Bis PAGE

Formulation and Storage

Formulation	Supplied as 0.22µm filtered solution in 50mM Tris, 150mM NaCl, 10% Glycerol (pH 8.0).
Storage	Valid for 12 months from date of receipt when stored at -80°C. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

D-dopachrome tautomerase (D-DT) is a newly described cytokine and a member of the macrophage migration inhibitory factor (MIF) protein superfamily. The D-DT protein also is present in most tissues and circulates in serum at similar concentrations as MIF. D-DT binds the MIF cell surface receptor complex, CD74/CD44, with high affinity and induces similar cell signaling and effector functions.

Assay Data

Tris-Bis PAGE



Cynomolgus DDT on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.