

Human DLL3 Domain (311-492) Protein, Ultra Low Endotoxin

Cat. No. DLL-HM1D4-UL

Description

Source	Recombinant Human DLL3 Domain (311-492) Protein is expressed from HEK293 with His tag at the C-terminus. It contains Val311-Leu492.
Accession	Q9NYJ7-1
Molecular Weight	The protein has a predicted MW of 20.26 kDa. Due to glycosylation, the protein migrates to 24-30 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE

Formulation and Storage


Formulation	Lyophilized from 0.22 µm filtered solution in PBS, 200mM L-Arginine (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Delta-like protein 3 (DLL3) is a transmembrane protein that belongs to the Delta/Serrate/Lag-2 (DSL) family of Notch ligands. DLL3 inhibits primary neurogenesis. May be required to divert neurons along a specific differentiation pathway. Plays a role in the formation of somite boundaries during segmentation of the paraxial mesoderm (By similarity).

Assay Data

Bis-Tris PAGE




Human DLL3 Domain (311-492) on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

ELISA Data

Assay Data

Human DLL3 Domain (311-492), His Tag ELISA


0.05µg Human DLL3 Domain (311-492), His Tag Per Well



Immobilized Human DLL3 Domain (311-492), His Tag at 0.5µg/ml (100µl/well) on the plate. Dose response curve for Anti-DLL3 Antibody, hFc Tag with the EC50 of 3.2ng/ml determined by ELISA.


Human DLL3 Domain (311-492), His Tag ELISA

0.05µg Human DLL3 Domain (311-492), His Tag Per Well



Immobilized Human DLL3 Domain (311-492), His Tag at 0.5µg/ml (100µl/well) on the plate. Dose response curve for Tarmatamab, hFc Tag with the EC50 of 2.8ng/ml determined by ELISA.

BLI Data



Loaded Anti-DLL3 Antibody, hFc Tag on ProA-Biosensor can bind Human DLL3 Domain (311-492), His Tag with an affinity constant of 0.62 nM as determined in BLI assay.