

Biotinylated Human ENPP-3 (48-157) Protein, Ultra Low Endotoxin



Cat. No. ENP-HM404B-UL

Description

Source	Recombinant Biotinylated Human ENPP-3(48-157) Protein is expressed from HEK293 with His tag and Avi tag at the N-Terminus. It contains Leu48-Asp157.
Accession	O14638
Molecular Weight	The protein has a predicted MW of 17.21 kDa. Due to glycosylation, the protein migrates to 18-20 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

Ectonucleotide pyrophosphatase-phosphodiesterase 3 (ENPP3), a protein detected in the human uterus, has been found to play an important role in the development and invasion of tumours. It was recently discovered that ENPP3 was upregulated during the window of implantation in the human endometrium but its functional relevance remains elusive.

Assay Data


Bis-Tris PAGE



Biotinylated Human ENPP-3 (48-157) on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Biotinylated Human ENPP-3 (48-157) is greater than 95% as determined by SEC-HPLC.