

Mouse EPHA2 Protein, Ultra Low Endotoxin



Cat. No. EPH-MM1A2-UL


Description	
Source	Recombinant Mouse EPHA2 Protein is expressed from HEK293 with His tag at the C-Terminus. It contains Lys26-Val538.
Accession	Q03145
Molecular Weight	The protein has a predicted MW of 57.95 kDa. Due to glycosylation, the protein migrates to 60-70 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background	
Erythropoietin-producing hepatocellular receptor A2 (EphA2) receptor tyrosine kinase plays an important role in tissue organization and homeostasis in normal organs. EphA2 is overexpressed in a variety of types of solid tumours with oncogenic functions.	

Assay Data


Bis-Tris PAGE



Mouse EPHA2 on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data



The purity of Mouse EPHA2 is greater than 95% as determined by SEC-HPLC.