

Mouse Fas/TNFRSF6/CD95 Protein

Cat. No. FAS-MM201

Description

Source	Recombinant Mouse Fas/TNFRSF6/CD95 Protein is expressed from HEK293 with hFc tag at the C-Terminus. It contains Gln22-Arg169.
Accession	P25446
Molecular Weight	The protein has a predicted MW of 43.3 kDa. Due to glycosylation, the protein migrates to 60-65 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 $\mu\text{g}/\text{ml}$ is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

CD95 (also known as Fas) is a member of the tumor necrosis factor receptor (TNFR) superfamily. Its cognate ligand, CD95L, is implicated in immune homeostasis and immune surveillance. Mutations in this receptor are associated with a loss of apoptotic signaling and have been detected in an autoimmune disorder called autoimmune lymphoproliferative syndrome (ALPS) type Ia, which shares some clinical features with systemic lupus erythematosus (SLE).

Assay Data


Bis-Tris PAGE



Mouse Fas on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Mouse Fas is greater than 95% as determined by SEC-HPLC.

ELISA Data

Mouse Fas, hFc Tag ELISA

0.1µg Human Fas Ligand, His Tag Per Well



Immobilized Human Fas Ligand, His Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Mouse Fas, hFc Tag with the EC50 of 0.25µg/ml determined by ELISA.