

# Cynomolgus/Rhesus macaque TNFSF15 Protein

Cat. No. FSF-CM115



## Description

<b>Source</b>	Recombinant Cynomolgus/Rhesus macaque TNFSF15 Protein is expressed from HEK293 with His tag at the N-Terminus.
	It contains Leu72-Leu251.
<b>Accession</b>	A0A2K5UA22
<b>Molecular Weight</b>	The protein has a predicted MW of 21.52 kDa. Due to glycosylation, the protein migrates to 25-35 kDa based on Tris-Bis PAGE result.
<b>Endotoxin</b>	Less than 1EU per µg by the LAL method.
<b>Purity</b>	> 90% as determined by Tris-Bis PAGE
	> 95% as determined by HPLC

## Formulation and Storage


<b>Formulation</b>	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
<b>Reconstitution</b>	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Storage</b>	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3-6 months after reconstitution. 2-8°C for 2-7 days after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

## Background

TNF superfamily member 15 (TNFSF15), a cytokine largely produced by vascular endothelial cells and a specific inhibitor of the proliferation of these same cells, can inhibit VEGF-induced vascular permeability in vitro and in vivo, and that death receptor 3 (DR3), a cell surface receptor of TNFSF15, mediates TNFSF15-induced dephosphorylation of VEGFR2.

## Assay Data


### Tris-Bis PAGE



Cynomolgus/Rhesus macaque TNFSF15 on Tris-Bis PAGE under reduced condition. The purity is greater than 90%.

### SEC-HPLC

Assay Data




The purity of Cynomolgus/Rhesus macaque TNFSF15 is greater than 95% as determined by SEC-HPLC.

Assay Data


ELISA Data

**Cynomolgus/Rhesus macaque TNFSF15, His Tag ELISA**  
0.2μg Cynomolgus/Rhesus macaque TNFSF15, His Tag Per Well



Immobilized Cynomolgus/Rhesus macaque TNFSF15, His Tag at 2μg/ml (100μl/well) on the plate. Dose response curve for Anti-TNFSF15 Antibody, hFc Tag with the EC50 of 3.0ng/ml determined by ELISA (QC Test).

SPR Data



Mouse DR3, His Tag immobilized on CM5 Chip can bind Cynomolgus/Rhesus macaque TNFSF15, His Tag with an affinity constant of 11.70 nM as determined in SPR assay (Biacore T200).