

Cat. No. FSF-HM418-UL

Description

Source	Recombinant Human GITR Ligand/TNFSF18 Trimer Protein is expressed from HEK293 with His tag and flag tag at the N-Terminus. It contains Gln50-Ser177 trimer design.
Accession	Q9UNG2
Molecular Weight	The protein has a predicted MW of 47.5 kDa. Due to glycosylation, the protein migrates to 54-64 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 0.01 EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC

Formulation and Storage

Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Dissolve the lyophilized protein in distilled water. Please refer to the Certificate of Analysis for detailed instructions.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background


Glucocorticoid-induced TNFR-related protein (TNFRSF18, GITR, CD357), expressed by T cells, and its ligand (TNFSF18, GITRL), expressed by myeloid populations, provide co-stimulatory signals that boost T cell activity. Due to the important role that GITR plays in regulating immune functions, agonistic stimulation of GITR is a promising therapeutic concept.

Assay Data**Bis-Tris PAGE**

Human GITR Ligand Trimer on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC

Assay Data




The purity of Human GITR Ligand Trimer is greater than 95% as determined by SEC-HPLC.

ELISA Data

Human GITR Ligand Trimer, His Tag ELISA

0.5 µg Human GITR Ligand Trimer, His Tag Per Well



Immobilized Human GITR Ligand Trimer, His Tag at 5 µg/ml (100 µl/well) on the plate. Dose response curve for Human GITR, hFc Tag with the EC50 of 50.8 ng/ml determined by ELISA.