# Human GARP (Y137H) & Latent TGF Beta 1 Complex Protein






Description	
Source	Recombinant Human GARP(Y137H)&Latent TGF Beta 1 Complex Protein is expressed from HEK293 with His tag and Avi tag at the C-Terminus.
	It contains His20-Leu628 (GARP(Y137H)) and Leu30-Ser390 (Latent TGF Beta 1).
Accession	Q14392(GARP(Y137H))&P01137(Latent TGF beta 1)
Molecular Weight	The protein has a predicted MW of 70.3 kDa (GARP(Y137H)) and 41.4 kDa (Latent TGF Beta 1). Due to glycosylation, the protein migrates to 71-75 kDa (GARP(Y137H)), 13 kDa and 42-45 kDa (Latent TGF Beta 1) based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	> 95% as determined by HPLC
Formulation and S	Storage
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt80°C for 3 months after reconstitution.Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

GARP&Latent TGF Beta is a complex found on surface of many types of cells. In Tregs, GARP is involved in TCR-mediated activation of Latent TGF- $\beta$  and thus promoting secretion and activation of TGF- $\beta$ . Integrin  $\alpha\nu\beta$ 8 on the surface of immune cells and other cells recognizes RGD in LAP, resulting in the release of mature TGF- $\beta$ 6 from the TGF- $\beta$ 8GARP complex.

## **Assay Data**

**Background** 

### **Bis-Tris PAGE**




Human GARP (Y137H) & Latent TGF Beta 1 Complex on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

**SEC-HPLC** 



### **Assay Data**




The purity of Human GARP (Y137H) &Latent TGF Beta 1 Complex is greater than 95% as determined by SEC-HPLC.

#### **ELISA Data**

### Human GARP (Y137H)&Latent TGF Beta 1, His Tag ELISA

0.5μg Human GARP (Y137H)&Latent TGF Beta 1, His Tag Per Well



Log Anti-GARP&TGF beta 1 Antibody, hFc Tag Conc.(μg/ml)

Immobilized Human GARP (Y137H) &Latent TGF Beta 1 Complex, His Tag at 5µg/ml (100µl/well) on the plate. Dose response curve for Anti-GARP&Latent TGF Beta 1 Antibody, hFc Tag with the EC50 of 14.8ng/ml determined by ELISA.