

Biotinylated Mouse GARP&Latent TGF beta 1 Complex Protein



Cat. No. GAT-MM401B

Description	
Source	Recombinant Biotinylated Mouse GARP&Latent TGF beta 1 Complex Protein is expressed from HEK293 with His tag and Avi tag at the N-Terminus. It contains Ala18-Leu628(GARP)&Leu30-Ser390(Latent TGF beta 1).
Accession	G3XA59(GARP)&P04202(Latent TGF beta 1)
Molecular Weight	The protein has a predicted MW of 66.6 kDa (GARP)&44.2 kDa (Latent TGF beta 1). Due to glycosylation, the protein migrates to 72-76 kDa (GARP)&13 kDa&45-47 kDa (Latent TGF beta 1) based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per µg by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC


Formulation and Storage	
Formulation	Lyophilized from 0.22µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

GARP&Latent TGF Beta is a complex found on surface of many types of cells. In Tregs, GARP is involved in TCR-mediated activation of Latent TGF-β and thus promoting secretion and activation of TGF-β. Integrin αβ8 on the surface of immune cells and other cells recognizes RGD in LAP, resulting in the release of mature TGF-β from the TGF-β&GARP complex.


Assay Data

Bis-Tris PAGE



Biotinylated Mouse GARP&Latent TGF beta 1 Complex on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.

SEC-HPLC



The purity of Biotinylated Mouse GARP&Latent TGF beta 1 Complex is greater than 95% as determined by SEC-HPLC.

Biotinylated Mouse GARP&Latent TGF beta 1 Complex Protein


Cat. No. GAT-MM401B

Assay Data

ELISA Data

Biotinylated Mouse GARP&Latent TGF beta 1, His Tag ELISA

0.1µg Mouse ITGAV&ITGB6, His Tag Per Well



Immobilized Mouse ITGAV&ITGB6, His Tag at 1µg/ml (100µl/well) on the plate. Dose response curve for Biotinylated Mouse GARP&Latent TGF beta 1 Complex, His Tag with the EC50 of 0.12µg/ml determined by ELISA.

Log Biotinylated Mouse GARP&Latent TGF beta 1, His Tag Conc.(µg/ml)