

Biotinylated Cynomolgus GIPR Protein

Cat. No. GIP-CM40RB

Description

Source	Recombinant Biotinylated Cynomolgus GIPR Protein is expressed from HEK293 with His tag and Avi tag at the C-terminus. It contains Gly26-Gln138.
Accession	XP_005589662.2
Molecular Weight	The protein has a predicted MW of 15.91 kDa. Due to glycosylation, the protein migrates to 30-40 kDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per μg by the LAL method.
Purity	>95% as determined by Bis-Tris PAGE >95% as determined by HPLC

Formulation and Storage


Formulation	Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 $\mu\text{g}/\text{ml}$ is recommended. Dissolve the lyophilized protein in distilled water.
Storage	-20 to -80°C for 12 months as supplied from date of receipt. -80°C for 3 months after reconstitution. Recommend to aliquot the protein into smaller quantities for optimal storage. Please minimize freeze-thaw cycles.

Background

The gastric inhibitory polypeptide receptor (GIPR), a G protein-coupled receptor (GPCR) that regulates glucose metabolism and insulin secretion, is a target for the development of therapeutic agents to address type 2 diabetes and obesity.

Assay Data

Bis-Tris PAGE



Biotinylated Cynomolgus GIPR on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.


SEC-HPLC

Biotinylated Cynomolgus GIPR Protein

Cat. No. GIP-CM40RB



Assay Data



The purity of Biotinylated Cynomolgus GIPR is greater than 95% as determined by SEC-HPLC.